



Excavator	42
Filling Instrument	44
Composite Syringe	47
Amalgam Product	47
Matrix Retainer / Separator	48
Contact Protector	48
Applicator	49




Excavator

-Round Tip Type-2-




	Description
Str.	Cavity preparation for wall and floor.
Cvd.	Used for undercut area, narrow rim of an opening and at the cavity angle.




-CDR Type-2-

*Excavator which pursued minimal intervention.



*Material / Si (Color Ring)








Restorative Treatment


Plastic Filling Instrument





Proximal Filling Instrument

*For efficient filling to proximal, thickness at the tip is thin.



Cement Filling Instrument




No. 03241



Composite Filling Instrument





^{*}Thickness / 0.2mm (#2, #3, #4, #5) 0.1mm (#6)

Composite Filling Instrument Tip (Replaceable Type)




*Remove all joint parts and replaceable tip before sterilization and disinfection.

*For more replaceable type products see P. 211.



-Handle-



-Tip (3 ea.)-









*Thickness / 0.2mm

Composite Syringe Type-2








No. 13251



Main Body

Nozzle & Plug S (Plug / Green)

No. 13252



*Resin of necessary quantity is filled in a nozzle, and a plug is inserted in and it is used by equipping into

No.	Description
13255	Main Body L
13256	Main Body S
13188	Nozzle & Plug L (100 ea.)
13190	Nozzle & Plug S (100 ea.)
13193	Nozzle & Plug L (500 ea.)
13194	Nozzle & Plug S (500 ea.)

*Material / PPS (Black part)

Amalgam Carver





No. 07100


Matrix Retainer Set








Individual Parts			
No.	Description		
13171	Main Body		
13172	Band (10 ea.)		



Separator





Contact Protector Set


Anterior (10 ea.) / Posterior (10 ea.)


*Stainless Steel

*Protect adjacent tooth during the cavity cutting.

No. 13175







Lining Applicator



-Round-





-Octagonal-





Ball-Pointed Applicator

No. 25015

