






Needle Holder	102
Gum Scissors	103
Tissue Pliers	104
Root Canal Instrument	106
Periosteal Elevator	108
Mirror	109









*Length is 180mm ,which is comfortable for microscope use.

Needle Holder MDS-180






Suture / from 6-0 to 8-0 Suture needle / from ø0.14 to ø0.33mm








Gum Scissors MDS-180







^{*}This feature can prevent the fingers from coming into sight.

Micro Tissue Pliers MDS-180







The precision of the tip



*Holding parts close horizontally and equally as surface-hold, not unequally like pint-hold. Therefore, you can surely hold delicate tissue and suture thread



*T/C is applied with tip part, which increases holding power of the Micro Tissue Pliers.




*Stopper, attached on the handle, always keeps the ideal distance. It prevents you from excessively opening and closing the pliers.

Tissue Pliers
























105

MTA Tweezers

toclavable
to 135°C

*For "Retrograde Filling of Root Canal".



Micro Root Canal Plugger





Micro Burnisher





Micro Explorer



*Effective in diagnosis in a bone cavity.



Micro Curette








Micro Periosteal Elevator






Micro Mirror Type-2



*Direct reflection type stainless mirror.

*Used with "Micro Handle"



3.0mm



*Even the part which is close to the edge is flat so those part also can reflect the object.

)	6.0mm	6.0mm	

2.0mm


ø5.0mm

No.	Description
14901	#1 (5 ea.)
14902	#2 (5 ea.)
14903	#3 (5 ea.)
14904	#4 (5 ea.)

[Clinical Movie]



Micro Handle



^{*}Easy detachable by twisting the end of the handle.



109



Micro Mirror Set Type-2



Contents Micro Handle / Micro Mirror Type2 #1, #2, #3, #4

No. 14900

