


# Elevator B Handle Spear

Suitably thin tip can approach a root and remove it from periodontal membrane sharply. This tip is able to hold a damage of surrounding tissues to a minimum so this elevator is effective for the extracting of "immediate implant placement".




Very sharp tip like a spear can find the point for approaching a root.




Very thin tip can approach the narrow space of periodontal membrane.

[How to use]




Put the tip into a root then periodontal membrane is amputated and the teeth is extracted smoothly.

[Attention]



Please do not lifting up a tooth by leverage-tequique. It may cause a tip broken.



# 1S  
No.16141



# 1C  
No.16144



# 2  
No.16145




# 3  
No.16146

※ maximum width of tip / 3.5 mm



[#1C]

The curve of the tip can approach the "molar" smoothly.



[#2 · #3]

This curved shape of a shank is suitably effective for "molar" to approach with ease.

The tip is positioned on the extended line of the handle for handling balance. It makes Dr,s comfortable to keep the balance of operation.

\*Because of YDM' s commitment to constant product improvement, YDM reserves the right to change prices and design or discontinue products without a notice.



MD516334 / ISO13485:2003  
FM516333 / ISO9001:2008  
Certification Site / Saitama Works

## YDM CORPORATION

6-5-20 Tabata, Kita-ku, Tokyo 114-0014, JAPAN

TEL:+81-3-3828-3161 FAX:+81-3-3827-8991

e-mail:ydm@ydm.co.jp

http://www.ydm.co.jp/en-new/index.html